

## Instructions for use of Isohelix SK2 Buccal Swabs







(REF SK-2S)

Intended for the retrieval of buccal cells. Single use only.



Use only if swab wrapper remains intact.

Follow the correct sampling procedure to avoid risk of the swab head becoming detached in the mouth. If the swab head does become detached from the swab stick whilst in the mouth, remove immediately.

<p><b>Take your DNA sample at least one hour after eating, drinking or cleaning your teeth.</b>  <b>For best results, rinse mouth with water immediately prior to sampling.</b></p>	
<p><b>1</b></p> 	<p>Pull open the package from one end.</p>
<p><b>2</b></p> 	<p>Remove one of the swabs from the tube.</p>
<p><b>3</b></p> 	<p>Insert the swab into your mouth and rub firmly against the inside of your cheek or underneath lower or upper lip. For standard DNA collection rub for <b>30 seconds</b> and in all cases rub for a minimum of 20 seconds.  <b>Important – use reasonable, firm and solid pressure</b></p>
<p><b>4</b></p> 	<p>Slide the plastic cap over the swab handle with the flat side of the cap facing upwards and the swab facing downwards.</p>
<p><b>5</b></p> 	<p>Insert the swab into the clear plastic tube and push the cap into place.</p> <p>Next, hold the cap while pulling the swab handle outwards to release the swab material into the tube.</p>
<p><b>6</b></p> 	<p>Close the cap.</p> <p>The tube is now completely sealed.</p>

May be used in conjunction with *In Vitro* Diagnostic Medical Devices. It is the responsibility of the Laboratory or Professional User to verify and validate the use of these DNA swabs with their device.

Isohelix is a division of Cell Projects

For Swab or DNA Isolation queries info@isohelix.com [www.isohelix.com](http://www.isohelix.com)

Molecular Biology Solutions [www.cellprojects.com](http://www.cellprojects.com)

Unit 2 Roebuck Business Park, Harrietsham, Kent, ME17 1AB, UK ☎: +44 1622 851177