

# Isohelix Swab and DNA Sample Handling Products Catalogue

Cat. No.	Description
<b>DNA Buccal Swabs</b>	
SK-1S	<b>100 x 1 swab</b> with 5ml tube and cap individually wrapped, and gamma treated.
SK-2S	<b>100 x 1 swab</b> with 2ml tube and SR- special release cap, individually wrapped and gamma treated.
SK-3S	<b>250 x 1 swab</b> individually wrapped and gamma treated.
SK-4S	<b>500 swabs</b> , 250 x 2 swabs per wrapper, gamma treated, packed 250 units to a box.
<b>DNA Isolation and Handling kits</b>	
DSK-3	1.5ml of Cell Lysis and Long term DNA stabilising buffer and PK tube for <b>3 reactions</b> .
DSK- 50	25ml of Cell Lysis and Long term DNA stabilising buffer and PK tube for <b>50 reactions</b>
DDK-3	Swab Optimised DNA Isolation kit for <b>3 reactions</b> containing 3 solutions.
<i>DDK-3/SK1</i>	As DDK-3 with 3 x SK-1 swabs
<i>DDK-3/SK2</i>	As DDK-3 with 3 x SK-2 swabs
DDK-50	Swab Optimised DNA Isolation kit for <b>50 reactions</b> containing 3 solutions.
<i>DDK-50/SK1</i>	As DDK-50 with 50 x SK-1 swabs
<i>DDK-50/SK2</i>	As DDK-50 with 50 x SK-2 swabs
BEK-3/SK1	Buccalyse DNA Extraction Kit – a quick and simple one tube method for extracting DNA from buccal swabs for 3 reactions. Supplied with 3 x SK-1 swabs
BEK-25/SK1	Buccalyse DNA Extraction Kit – a quick and simple one tube method for extracting DNA from buccal swabs for 25 reactions. Supplied with 25 x SK-1 swabs
DQC-10	Pre check for quality test of Human DNA taken from swabs for <b>10 reactions</b>
DQC-50	Pre check for quality test of Human DNA taken from swabs for <b>50 reactions</b>
SGC-3/SK1	<b>3x</b> Isohelix Dri-Capsules for long term DNA stability on SK-1 swabs. Individually foil wrapped.
SGC-50/SK1	<b>50x</b> Isohelix Dri-Capsules for long term DNA stability on SK-1 swabs. Individually foil wrapped.